

## NEW VISUALIZATION SYSTEMS FROM FUJIFILM



**FUJIFILM New Development, USA**, has introduced 2 visualization systems for minimally invasive surgery. Using proprietary technology, the **Ultra-Slim Video Laparoscope System (EL-580FN)** delivers enhanced image resolution, color fidelity, and display quality, says **FUJIFILM**. The product features include “Chip on the Tip” high-definition digital imaging

processing, less fogging, autoclave sterilization reprocessing, and a low profile, lightweight ergonomic handle. The 3.8-mm-diameter distal end was designed to improve workflow, reduce physician fatigue, and potentially reduce the size of incisions. The accompanying Digital Video Processor System is used for endoscopic procedures with automatic light control, an anti-blur function for motion images, and digital zoom.

**FUJIFILM** reports that the **Full High Definition Surgical Visualization System** is designed for a wide variety of surgical applications and offers edge enhancement, automatic gain control, dynamic contrast function, selective color enhancement, smoke reduction, and grid removal features. It includes a portfolio of rigid scopes, cameras, and video processing systems.

**FOR MORE INFORMATION, VISIT:** <http://www.fujifilmusa.com>

## FREEMIE BREAST MILK COLLECTION SYSTEM



**Freemie**® offers a hands-free breast-milk collection system with the **Freemie Liberty Mobile Hands Free Breast Pump**

**System** and **Next Generation Freemie Closed System Collection Cups**.

The concealable pump has a rechargeable battery and hospital-power suction for single or double pumping. Programmable memory buttons allow the mother to pre-set or adjust speed and suction functions. Tubing lengths can be changed so that the pump can be placed on a desk, worn with a detachable belt clip, or carried in a bag.

**Freemie** says the cups are lower-profile and more compact than other pump system cups, and when placed on the breast under the mother’s bra, can be easily removed so that milk can be transferred to storage. Each cup, with a 25 mm or 28 mm funnel and valve, holds 8 oz of milk.

**FOR MORE INFORMATION, VISIT:** <http://www.freemie.com>

## PREVENTEZA: EMERGENCY CONTRACEPTIVE



**Combe, Inc**, the maker of **Vagisil**®, has launched **Preventeza™** (levonorgestrel tablet, 1.5 mg), an emergency contraceptive for the prevention of pregnancy if unprotected sex or failed birth control occurs.

Available online or over-the-counter as a single tablet, **Preventeza** is a proven option to help women prevent pregnancy

before it starts by using a higher dose of levonorgestrel than most birth control pills. It must be used within 72 hours of unprotected intercourse, and is not intended to be used as regular birth control. **Combe** says that **Preventeza** works mainly by stopping the release of an egg from the ovary and may also prevent fertilization of an egg or prevent a fertilized egg from implanting in the uterus. **Combe** also says that levonorgestrel 1.5 mg will not work if the woman is already pregnant and will not affect an existing pregnancy.

**FOR MORE INFORMATION, VISIT:** <https://www.vagisil.com/products/preventeza-emergency-contraceptive>

## UPSPRING'S C-PANTY FOR POSTCESAREAN RECOVERY



**UpSpring**® says that its patented **C-Panty**® undergarment provides medical-grade compression and speeds recovery after cesarean delivery. **C-Panty** helps to reduce swelling and discomfort, supports weakened muscles, and reduces the incision bulge without hooks, straps, or Velcro that might irritate the incision area.

**C-Panty's** medical-grade silicone panel suppresses the formation of excess or improperly formed collagen, which can contribute to scarring, says **UpSpring**. The silicone may help reduce itchiness, and can lessen the chance of infection at the incision area. The silicone is durable, washable, and integrated into the panty, eliminating the need for scar gel or scar gel pads.

The **C-Panty** can be worn immediately after birth and for up to 12 months. If worn when the incision is not healed, the silicone panel should be covered with a panty liner or pad. Once the incision has healed, the covering can be discontinued.

**FOR MORE INFORMATION, VISIT:** <https://www.upspringbaby.com/cpanty>